





---

## Intra-Aortic Balloon Pump

### MCare Service Plans

---

**In addition to the standard warranty provided with your IABP**, Maquet offers three different levels of Extended Warranty Programs and Service Contracts for additional protection.

The Extended Warranty program is available at the time of the initial purchase or within the first six months of delivery of your IABP.

If a service need arises within the initial 12 month standard warranty period, Maquet will provide repair services as necessary at no charge. Maquet warrants that components within our products will be free from all defects in workmanship and materials.

Service Plan options are available for any previously purchased products without extended warranty coverage.



## Cables & Leadwires

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0012-00-1155-01	ECG Patient Cable - 5 Lead - 10' (304.8cm) AAMI	✓	✓
0012-00-1155-02	ECG Patient Cable - 5 Lead - 10' (304.8cm) IEC	✓	✓
0012-00-1156-01	ECG Leadwires - 5 Lead - Snap Type - 50" (127cm) AAMI	✓	✓
0012-00-1156-02	ECG Leadwires - 5 Lead - Snap Type - 50" (127cm) IEC	✓	✓
0012-00-1157-01	ECG ESIS Leadwires - 5 Lead - Pinch Type - 50" (127cm) AAMI	✓	✓
0012-00-1157-02	ECG ESIS Leadwires - 5 Lead - Pinch Type - 50" (127cm) IEC	✓	✓
0012-00-1212	Remote Monitor Extension Cable - 10' (304.8cm)	✓	✓
0012-00-1245	Interface Transducer Cable - 12' (365.7cm)	✓	✓
0012-00-1589-01	Low Level Output Cable - Underminated - 15' (457.2cm)	✓	
0012-00-1589-02	Low Level Output Cable - Interface to Datascope Monitor - 15' (457.2cm)	✓	
0012-00-1589-03	Low Level Output Cable - Interface to Philips monitor - 14' (426.7cm)	✓	
0012-00-1589-04	Low Level Output Cable - Interface to GE monitor - 14' (426.7cm)	✓	
0012-00-1589-05	Low Level Output Cable - Underminated - 23' (701cm)	✓	



ECG Patient Cable - 5 Lead - 10' (304.8cm) - AAMI  
0012-00-1155-01



ECG Patient Cable - 5 Lead - 10' (304.8cm) - IEC  
0012-00-1155-02



ECG Leadwires - 5 Lead - Snap Type - 50" (127cm) - AAMI  
0012-00-1156-01



ECG Leadwires - 5 Lead - Snap Type - 50" (127cm) - IEC  
0012-00-1156-02



ECG ESIS Leadwires - 5 Lead - Pinch Type - 50"  
(127cm) - AAMI 0012-00-1157-01



ECG ESIS Leadwires - 5 Lead - Pinch Type - 50"  
(127cm) - IEC 0012-00-1157-02



Remote Monitor Extension Cable - 10' (304.8cm)  
0012-00-1212



Interface Transducer Cable - 12' (365.7cm)  
0012-00-1245



Low Level Output Cable - Unterminated  
15' (457.2cm) 0012-00-1589-01  
23' (701cm) 0012-00-1589-05



Low Level Output Cable - Interface to  
Datascope Monitor - 15' (457.2cm)  
0012-00-1589-02



Low Level Output Cable - Interface to  
Philips Monitor - 14' (426.7cm) 0012-00-1589-03



Low Level Output Cable - Interface to  
GE Monitor - 14' (426.7cm) 0012-00-1589-04

## Cables & Leadwires (continued)

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0012-00-0323	External Signal Cable - Unterminated - 15' (457.2cm)	✓	✓
0012-00-0819-01	Cable Interfaces BP & ECG for HP/Philips Monitors, except CMS, Viridia or Intellivue - 20' (609.6cm)	✓	✓
0012-00-0819-02	Cable Interfaces BP & ECG for HP/Philips Monitors, except CMS, Viridia or Intellivue - 30' (914.4cm)	✓	✓
0012-00-0819-03	Cable Interfaces BP & ECG for HP/Philips Monitors, except CMS, Viridia or Intellivue, RoHS - 30' (914.4cm)	✓	✓
0012-00-0820-03	Cable Interfaces ECG for HP/Philips CMS, Viridia, Intellivue & MP series monitors, RoHS - 30' (914.4cm)	✓	✓
0012-00-0821-02	Cable Interfaces BP & ECG for Siemens monitors with 1/4" stereo phone plugs - 30' (914.4cm)	✓	✓
0012-00-0822-02	Cable Interfaces BP for Philips Viridia & BP & ECG for Spacelabs PC series monitors - 30' (914.4cm)	✓	✓
0012-00-0822-03	Cable Interfaces BP for Philips Viridia & BP & ECG for Spacelabs PC series monitors, RoHS - 30' (914.4cm)	✓	✓
0012-00-0823-02	Cable Interfaces BP for Marquette 7000 series monitors - 30' (914.4cm)	✓	✓
0012-00-0824-02	Cable Interfaces BP & ECG for Marquette 7000 series monitors - 30' (914.4cm)	✓	✓
0012-00-0825-02	Cable Interfaces BP & ECG for GE Marquette monitors - 30' (914.4cm)	✓	✓
0012-00-0825-03	Cable Interfaces BP & ECG for GE Marquette monitors, RoHS - 30' (914.4cm)	✓	✓
0012-00-0826-03	Cable Interfaces ECG for Spacelabs PC series monitors with TT-253 plugs, RoHS - 30' (914.4cm)	✓	✓
0012-00-1650-01	Cable Interfaces Datascope Spectrum/OR monitors to ECG & IBP Sync of IABP - 15' (457.2cm)	✓	✓



External Signal Cable - Unterminated  
15' (457.2cm) 0012-00-0323



Cable Interfaces BP & ECG for HP/Philips Monitors, except CMS, Viridia or Intellivue - 20' (609.6cm)  
0012-00-0819-01



Cable Interfaces ECG for HP/Philips CMS, Viridia, Intellivue & MP series monitors , RoHS - 30' (914.4cm) 0012-00-0820-03



Cable Interfaces BP & ECG for Siemens monitors with 1/4" stereo phone plugs - 30' (914.4cm) 0012-00-0821-02



Cable Interfaces BP for Philips Viridia & BP & ECG for Spacelabs PC series monitors - 30' (914.4cm) 0012-00-0822-02



Cable Interfaces BP for Marquette 7000 series monitors - 30' (914.4cm) 0012-00-0823-02



Cable Interfaces BP & ECG for Marquette 7000 series monitors - 30' (914.4cm) 0012-00-0824-02



Cable Interfaces BP & ECG for GE Marquette monitors - 30' (914.4cm) 0012-00-0825-02



Cable Interfaces ECG for Spacelabs PC series monitors with TT-253 plugs, RoHS - 30' (914.4cm) 0012-00-0826-03



Cable Interfaces Datascope Spectrum/OR monitors to ECG & IBP Sync of IABP - 15' (457.2cm) 0012-00-1650-01

## AC Line Cords

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0012-00-0886-01	AC Line Cord - 110v - United States - 16' (4.9m)	✓	✓
0012-00-0886-02	AC Line Cord - 220v - International - 16' (4.9m)	✓	✓
0012-00-0886-03	AC Line Cord - 220v - China - 16' (4.9m)	✓	✓
0012-00-0886-04	AC Line Cord - 220v - Brazil - 16' (4.9m)	✓	✓



AC Line Cord - 110v - United States - 16' (4.9m)  
0012-00-0886-01



AC Line Cord - 220v - International - 16' (4.9m)  
0012-00-0886-02



AC Line Cord - 220v - China - 16' (4.9m)  
0012-00-0886-03



AC Line Cord - 220v - Brazil - 16' (4.9m)  
0012-00-0886-04



## Safety Disk

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0997-00-0985-01	Safety Disk Assembly - English	✓	✓
0997-00-0985-03	Safety Disk Assembly - German	✓	✓
0997-00-0985-04	Safety Disk Assembly - French	✓	✓
0997-00-0985-05	Safety Disk Assembly - Spanish	✓	✓
0997-00-0985-06	Safety Disk Assembly - Japanese	✓	✓
0997-00-0985-07	Safety Disk Assembly - Italian	✓	✓
0997-00-0985-08	Safety Disk Assembly - Portuguese	✓	✓
0997-00-0985-09	Safety Disk Assembly - Swedish	✓	✓
0997-00-0985-11	Safety Disk Assembly - Danish	✓	✓
0997-00-0985-15	Safety Disk Assembly - Chinese	✓	✓



Safety Disk Assembly 0997-00-0985-01

## Condensate Removal Module

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0997-00-0986-01	Condensate Removal Module - English	✓	✓
0997-00-0986-03	Condensate Removal Module - German	✓	✓
0997-00-0986-04	Condensate Removal Module - French	✓	✓
0997-00-0986-05	Condensate Removal Module - Spanish	✓	✓
0997-00-0986-06	Condensate Removal Module - Japanese	✓	✓
0997-00-0986-07	Condensate Removal Module - Italian	✓	✓
0997-00-0986-08	Condensate Removal Module - Portuguese	✓	✓
0997-00-0986-09	Condensate Removal Module - Swedish	✓	✓
0997-00-0986-11	Condensate Removal Module - Danish	✓	✓
0997-00-0986-15	Condensate Removal Module - Chinese	✓	✓



Condensate Removal Module 0997-00-0986-01

## Helium Tanks & Parts

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0075-02-0001-03	Helium Tank - UK / France - Reusable - Pack of 3 Tanks	✓	✓
0075-02-0002-03	Helium Tank - German - Reusable - Pack of 3 Tanks	✓	✓
0075-00-0024-01	Helium Tank - Reusable - Single Tank	✓	✓
0075-00-0024-03	Helium Tank - Reusable - Pack of 3 Tanks	✓	✓
0075-00-0034-03	Helium Tank - Reusable - Pack of 3 Tanks - Empty	✓	✓
0202-00-0104	Helium Tank - Disposable - Single Tank - Full	✓	✓
0348-00-0185	Washer - Helium Tank, 1 pc.	✓	✓
0366-00-0092	Valve Knob - Helium Tank - United States	✓	✓
0366-00-0098	Valve Knob - Helium Tank - Metric	✓	✓



Helium Tank - UK / France - Reusable - Pack of 3 Tanks 0075-02-0001-03



Helium Tank - German - Reusable - Pack of 3 Tanks 0075-02-0002-03



Helium Tank - Reusable - Single Tank 0075-00-0024-01



Helium Tank - Reusable - Pack of 3 Tanks - Empty 0075-00-0034-03



Helium Tank - Disposable - Single Tank - Full 0202-00-0104



Washer - Helium Tank, 1 pc. 0348-00-0185



Valve Knob - Helium Tank - United States 0366-00-0092

## Accessories & Replacement Parts

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0020-00-0418-01	Manual Fill Kit - English	✓	✓
0146-00-0039	12v Battery - Single Battery Only - Sealed Lead Acid - Replace in Pairs	✓	✓
0146-00-0047-01	Battery - Complete Assembly - Sealed Lead Acid	✓	✓
0154-01-0001	Doppler - Complete Assembly	✓	✓
0154-01-0005	Doppler - Probe Only	✓	✓
0198-00-0023	Protective Cover		✓
0380-00-0284	Guard - Condensate Remove Module	✓	✓
0406-00-0742	Replacement IV Pole Mounting Bracket	✓	✓
0436-00-0199	Replacement IV Pole - Pole Only	✓	✓
0997-00-0417-01	Remote Monitor Mount	✓	✓
0997-00-0428	Assembly, Accessory Storage Pouch		✓
0997-01-0532	Storage Bag (without Doppler)	✓	



Manual Fill Kit - English  
0020-00-0418-01



12v Battery - Single Battery  
Only - Sealed Lead Acid  
0146-00-0039  
Replace in pairs



Battery - Complete Assembly  
Sealed Lead Acid  
0146-00-0047-01



Doppler - Complete Assembly  
0154-01-0001




Doppler - Probe Only  
0154-01-0005



Protective Cover 0198-00-0023



Guard - Condensate Remove  
Module 0380-00-0284



Replacement IV Pole Mounting  
Bracket 0406-00-0742



Replacement IV Pole - Pole Only  
0436-00-0199



Remote Monitor Mount  
0997-00-0417-01



Assembly, Accessory Storage  
Pouch 0997-00-0428

## Consumable & Disposable Items

Part No.	Catalog Description	Compatability	
		CS300	CS100
0103-00-0167	3 Way Stopcock	✓	✓
0683-00-0422-02	Recorder Paper - box of 10 rolls	✓	✓



3 Way Stopcock 0103-00-0167



Recorder Paper 0683-00-0422-02  
Package of 10 rolls

## Maintenance Parts & Special Tools

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0040-00-0146	Preventive Maintenance Kit - 2,500 Hour	✓	✓
0040-00-0147	Preventive Maintenance Kit - 5,000 Hour	✓	✓
0040-00-0437	Replacement FOS Lamp Kit	✓	
0453-00-0154	Calibrated Syringe, 38.5cc	✓	✓
0683-00-0314	Calibration Chamber	✓	✓
0992-00-0245	FOS Interface Test Chamber	✓	
KIT-CLEANING	Fiber-Optic Cleaning Kit	✓	



Preventive Maintenance Kit - 2,500 Hour  
0040-00-0146



Preventive Maintenance Kit - 5,000 Hour  
0040-00-0147



Replacement FOS Lamp Kit 0040-00-0437



KIT-CLEANING Fiber-Optic Cleaning Kit

## Service Manuals

Part No.	Catalog Description	Compatibility	
		CS300	CS100
0070-00-0528-01	Service Manual - CS100		✓
0070-CD-0528-01	CD Service Manual - CS100		✓
0070-00-0689	Service Manual - CS300	✓	
0070-CD-0689	CD Service Manual - CS300	✓	











---

**MAQUET**  
GETINGE GROUP

Maquet  
1300 MacArthur Blvd., Mahwah, NJ 07430, USA  
Phone: +1 201 995 8700 or 1 800 777 4222  
Fax: +1 201 995 8700  
Email: [CA.Techsupport@maquet.com](mailto:CA.Techsupport@maquet.com)

[www.maquet.com](http://www.maquet.com)  
<http://ca.maquet.com>

---

**GETINGE GROUP**

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, Getinge and Maquet. ArjoHuntleigh focuses on patient mobility and wound management solutions. Getinge provides solutions for infection control within healthcare and contamination prevention within life sciences. Maquet specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.